
Week 2: Oscilloscopes and Function Generators

RMS of a periodic signal is calculated by first squaring the waveform, then taking its mean over its period, T, then taking the square root. Its definition using the calculus is

RMS =
$$\sqrt{\frac{1}{T} \int_0^T f^2(t) dt}$$

As an example, we will derive the equation for RMS/Vpp for a sine wave. You will be asked to derive the equation for square waves and triangular waves in the pre-lab.



First,
$$f(t) = a \sin\left(\frac{2\pi t}{T}\right) \to RMS = \sqrt{\frac{1}{T} \int_0^T a^2 \sin^2\left(\frac{2\pi t}{T}\right) dt}$$

Using the definition of
$$\sin^2(\theta) = \frac{1 - \cos(2\theta)}{2}$$
, $RMS = \sqrt{\frac{a^2}{T} \int_0^T \frac{1 - \cos(\frac{4\pi t}{T})}{2} dt}$

Taking the integral,
$$RMS = \sqrt{\frac{a^2}{T} \left[\frac{1}{2} t - \frac{\frac{1}{2}T}{4\pi} \sin\left(\frac{4\pi t}{T}\right) \right]_0^T}$$

Evaluating, we get
$$RMS = \sqrt{\frac{a^2}{T} \left[\frac{1}{2} T \right]}$$
 (Note that at $t = 0$ and T, $\sin \left(\frac{4\pi t}{T} \right) = 0$)

Therefore,
$$RMS = \frac{a}{\sqrt{2}}$$
, and since $Vpp = 2a$, then $\frac{RMS}{Vpp} = \frac{1}{2\sqrt{2}}$

It may for the purposes of your lab helpful to think of RMS in terms of Vpp, like so:


$$RMS = \frac{Vpp}{2\sqrt{2}}$$

Week 2 Prelab

Calculate the ratio RMS/Vpp for the following signals. Show all your work! Name:


1. Square Wave: RMS / Vpp = ?

UID: 705575353



Vpp = 2a $RMS = \int \int \int \int f^{2}(t) dt$ $= \int \int \int d^{2}t dt$ = 0

2. Triangular Wave: RMS / Vpp = ?



VPP = 2a

PMS: \[\frac{1}{7} \int \text{T'(1) det} \]

= \frac{1}{7} \fra

3. If you see a difference by a factor of 10 between the oscilloscope reading and the function generator setting, where is the first place that you should look? Watch the Probe Setting video (https://youtu.be/dtSuTHIviSo) for the answer.

the first place to look is the probe attenuation and the channel ratio.

4.	If you see a difference by a factor of 2 between the oscilloscope reading and the
	function generator setting, where is the first place that you should look? Watch the
	Function Generator Output Impedance video (https://youtu.be/-8Dv1oOjD9w) for
	the answer.

Look at the utility botton -> channel 1. There might be intenal pesistance.

5. Why would you ever want to use AC coupling on an oscilloscope? Watch the AC Coupling video (https://www.youtube.com/watch?v=dtSuTHIviSo&t=6s) for the answer.

Let you incresse the sestivity and see noise and then leep moresing until you see structure.

Week 2 Prelab End

Time Dependent Measurements

This week's experiments will give you the opportunity to learn the basic operations of an oscilloscope and a function generator.

Setting up Function Generator and Oscilloscope

- 1. Turn on both the function generator (Figure 2-1a) and oscilloscope (Figure 2-2).
- 2. Connect the function generator's CH1 output to the CH1 input of the oscilloscope.
- 3. Press the function generator Output button next to the CH1 connector if it is not lit.
- 4. You will now learn to display the input signal properly by using the three basic functions of the oscilloscope: Horizontal control, Vertical control, and Triggering.

5. [INITIALIZATION]

- a. Press the Save/Recall button in the Measure section of the controls.
- b. Push the Recall soft switch.
- c. Push the "Load from setup_0" soft switch.
- d. Press the "Press to Recall" soft switch. This is your starting point.
- 6. The display should be similar to that shown in Figure 2-1.




FIGURE 2-1. .Unstabilized, Horizontally and Vertically Incorrect, Sinusoid Display.